

## Hardware Field Guide

# Nebula Capsule Max



## Field Guide Contents

- 3 Your Invincikids kit
- 3 BERT setup
- 4 Projector controls
- 4 Powering on
- 5 Sound
- 5 Content
- 6 Handling and care

# Commitment to safety

Most importantly, we value your and your patients' safety and well-being. As with all media devices and new therapy tools, first familiarize yourself with the benefits, risks, and best practices.

## About this guide

### Technical user manual

This guide is not a replacement for the official device user manual. Refer to the manufacturer's manual for all headset safety guidelines and technical instructions.

**Refer the Nebula Capsule Max User Manual for technical instructions.**

[!\[\]\(3211b5d1d968fc1665909b34f9f16010\_img.jpg\) Nebula Capsule Max Manual](#)

The Nebula Capsule Max is manufactured by Anker.

### Annotations

Throughout the field guide, we've sprinkled in tips to help you optimize your projector use and patient experience.

Lookout for this  symbol for helpful user tips and tricks!



# Your kit

---

Your Nebula Capsule projector is part of the Invincikids **Bedside Entertainment and Relaxation Theater (BERT)** along with a screen and mounting hardware.

The BERT offers a way to entertain and distract bed-bound patients during potentially stressful moments. It includes experiences specifically designed to help patients through anesthesia inductions and for in-bed entertainment.

## Your kit includes:

1. Nebula Capsule Max projector
2. Projector mount
3. Screen mount
4. Nebula remote




# Device

---

## BERT Setup

The BERT requires some setup before use. First, mount the projector to the head of the patient's bed or to a gurney rail.


Then, mount the screen at the foot of the bed using the provided clamps or slide the mount under the mattress and elevate the screen, depending on your model.



# Projector controls


The top of the projector is a touch-sensitive control pad. This can be used to:

- Select apps
- Control induction experiences
- Adjust volume
- Play, pause, or fast-forward movies
- Exit apps




## The remote

The remote signal is detected by the sensor on the back of the projector. To use the remote like a computer mouse, press the button on the remote with the mouse icon. Move the cursor using the arrows on the remote and the OK button to select.



# Powering on

Press and hold the power button on the rear of the projector.



**To select:** Use the touch controls on the top of the projector or the OK button on the remote.



**Do not worry about focusing the projector—it's automatic.** If you want to manually focus, press and hold the option button on the remote (≡) for three seconds. Use the left and right arrows to adjust the focus, then press OK.

# Sound

The body of the projector is a built-in speaker. The projector can be paired with bluetooth headphones or with a smart device to function as an external speaker.

# Content

## Invincikids Induction Experiences

Your Nebula projector comes preloaded with experiences designed by Invincikids specifically to make mask inductions less scary and more fun.



Example Invincikids experience: Sevo the Dragon



Example Invincikids experience: Sevo & Desi



Invincikids Induction Experiences work particularly well with patients ages 3-8 who are engaged and willing to connect with the staff.

These are interactive experiences, guided by the provider. They require an enthusiastic provider who can craft an experience that meets a patient's individual needs.

**The provider must be present the entire time to control the game.**

## Available Commercial Apps and Services

### Movies Anywhere

The Nebula projector comes preloaded with family-friendly movies that are accessible offline, with the option to purchase additional movies through Movies Anywhere.

### Streaming services

Streaming services are available when the Nebula projector is connected to the internet via WiFi. These can be added through the Nebula Manager and require a login to access content.

### Third-party games

You can download third-party games from the Nebula Manager app on the projector home screen. Commercial apps will require your own account and login.



Streaming content requires a WiFi connection. If you are moving a patient between rooms, like traveling from pre-op to the operating room, ensure you have overlapping network coverage to avoid a lapse in video playback.



# Handling and care

---

## Cleaning

After every use, clean the projector and remote with disinfecting wipes. Avoid spraying liquid cleaners directly on the projector.